

Wallaby®3 Phototherapy System

Ideal for hospital and home use.



Children's
Medical
Ventures

RESPIRONICS®

Wallaby®3 phototherapy system allows a mother to hold, even nurse, her baby without interrupting therapy. In fact, at just 7.2 pounds, it's light enough that it can be easily carried around with the baby. That's a much more natural way to deliver the precise level of therapy. And it gives new mothers more time to bond with their babies.

Wallaby[®]3 Phototherapy System

The Wallaby 3 phototherapy system is a lightweight and compact alternative for delivering effective therapy to jaundiced babies. Every aspect of this device has been designed to make it the most complete – and most portable – phototherapy system available.

Efficient

FIBER OPTIC TECHNOLOGY:

- Softer, brighter pattern of therapeutic light – increased irradiance improves therapy delivery.

Safe

HOUR METER:

- Automatically records treatment time.

ENHANCED SAFETY FEATURES:

- Locking device prevents the panel from disengaging from the illuminator if the cable is pulled.
- Light shuts off automatically if fiber optic cable is disconnected.
- No infrared or ultraviolet light; not harmful to baby's eyes.

Flexible

TWO IRRADIANCE LEVELS:

- Increased light intensity to reduce higher bilirubin levels faster.

TWO-LAMP CAROUSEL: with Automatic Bulb Exchange

- Automatically exchanges the expired bulb, eliminating disruption of therapy and after-hours calls. A visual indicator alerts the caregiver that the second bulb is in use.

TWO SIZES OF FIBER PANELS:

- Choose the optimal modality for the infant: wrapped around the baby's torso or placed underneath the baby's body.



This neonatal phototherapy panel and the larger 14-inch wrap-around panel give you options to better tailor therapy to an individual patient's needs.

ordering info

item no.	description	quantity
H-3600	Wallaby 3 Phototherapy System	1 each
CC302	Wallaby 3 Carrying Case	1 each
EG-2000	Wrap-Around Panel	1 each
EG-2000N	NeoNatal Panel	1 each
JD-101	Joey Dosimeter	1 each
MD-510	Disposable Covers For Wrap-Around Panel	50/pkg
MD-610	Disposable Covers and Vest For NeoNatal Panel	100/pkg
PGM302	Parents' Guide	1 each
HPK-301	Patient Kit	1 each
H100-3600-20	Replacement Phototherapy Light Bulb	1 each
161-5000-00	Replacement Power Cord	1 each

specifications

ILLUMINATOR SIZE	0.25" x 9" x 4" (369 in.3); (26 cm x 23 cm x 10 cm (5980cm ³))	
WEIGHT	7.2 lbs (3.3 kg)	
FIBER OPTIC PANELS		
NeoNatal Panel	4" x 5" (10 cm x 13 cm)	
Standard Wrap-Around Panel	3" x 14" (7.5 cm x 35.5 cm)	
IRRADIANCE LEVELS	LEVEL 1	LEVEL 2
NeoNatal Panel	25-35 μ W/cm ² /nm	50-60 μ W/cm ² /nm
Standard Wrap-Around Panel	10-20 μ W/cm ² /nm	15-25 μ W/cm ² /nm



Customer Service: 1-800-345-6443 or 724-387-4000

Respironics Europe: +33-1-55-60-19-80

Respironics Asia Pacific: +81-3-5280-9611

Wallaby is a registered trademark of Respironics, Inc. People. Products. Programs. is a servicemark and Power Programs is a trademark of Respironics, Inc. and its affiliates. ©2004 Respironics, Inc.

CAUTION: U.S. federal law restricts these devices to sale by or on the order of a physician.

1024141 KB 10/4/04 G

